

CI 21

NH₃ Transmitter

Operations Manual



GfG Instrumentation

1194 Oak Valley Dr, Ste 20, Ann Arbor MI 48108 USA (800) 959-0329 • (734) 769-0573 • www.goodforgas.com/us-en

Table of Contents


	Page
For Your Safety	3
General Description	3
Detection Principle	4
Operation	4
Design	4
Mounting	5
Electrical Connections	6
Start-up	6
Check Electrical Zero Point	7
Check Sensitivity (Calibration) / Bump Test	7
Check Zero Point (Dead Band)	8
Troubleshooting	8
Service	8
Maintenance and Inspection	9
Spare Parts	9
Accessories	9
Connection Diagram Cl 21 (Internal Sensor)	10
Connection Diagram Cl 21 (External Sensor)	
PLC Charts	
Technical Data	14
Warranty	14

For Your Safety

As with any piece of complex equipment, the GfG CI 21 transmitter will do the job it is designed to do only if it is used and serviced in accordance with the manufacturer's instructions. Please protect yourself and your employees by following the instructions in this manual. All individuals who have or will have the responsibility for using and servicing this product must carefully read this manual. The warranties made by GfG with respect to the product are voided if functions or parameters are changed without the permission of GfG. They are also voided if the product is not used and serviced in accordance with the instructions in this manual. Failures or false alarms caused by interfering gases or electrical signals are not part of the warranty obligation. The above does not alter any statements by GfG regarding warranties, conditions of sale and/or delivery.

General Description

A fixed gas monitoring system consists of the transmitter (MWG) and a controller (GMA), which are connected by a cable. The transmitter converts the unit of the gas concentration into an electrical signal and transmits it to the controller.



Fixed Gas Detection System

Figure 1

The CI 21 is an "intelligent" transmitter. The comprehensive electronics allow easy operation and maintenance and also increase operational safety and accuracy. The following features characterize the CI 21:

- Easy on-site calibration
- Compensation for temperature fluctuations
- Long sensor life

Detection Principle


The CI 21 has a Charge Carrier Injection sensor. As long as there is no gas in the ambient air, the internal resistance of the sensor is high. NH₃ (ammonia) is adsorbed at the sensor surface, reducing the internal resistance. The change in resistance is equivalent to the measurement of the gas concentration and is converted into a standardized current signal (0.2–1 mA or 4–20 mA) by the integrated electronics. The Charge Carrier Injection sensor has proven long-term stability.

Operation

The CI 21 passes function and performance tests before shipment and is calibrated with suitable test gases. This does not preempt the user's obligation to perform a function test after putting the transmitter into operation.

Design

The design of the CI 21 is shown in figure 2. The sensor is mounted in the sensor support (pos. 5). Gas enters the sensor chamber through the diffusion inlet (pos. 6). The electronics convert the measurement signal into the output of 0.2–1 mA or 4–20 mA, which is then transferred to the controller by a cable.



4




Figure 3 / Schematic

Mounting

It is essential to know the ambient conditions before deciding on the mounting position. To achieve representative measurement results, check the:

- Room ventilation
- Gas density

Install the CI 21 in a place where gases will pass the sensor even in case of bad ventilation. If necessary, use a smoke cartridge to check air flow.

Take the following into consideration as well:

- Ambient temperature
- Rain water, hose water, dripping water or condensation
- Dust in the atmosphere
- Avoid installation near silicones, anti-seize, lubricants
- Other possible contaminants

The CI 21 is protected to a great extent against water and dust (IP 54). Special accessories are available to provide additional protection for very harsh conditions. Please contact GfG for detailed information.



The warranty may be voided if the transmitter is exposed to ambient conditions that were unknown to GfG during planning, production or delivery.

Substitution of components may compromise the safety of this device. Only use components approved by GfG Instrumentation, Inc.

When deciding on the transmitter position, make sure that it is always accessible for service and maintenance. Mounting the CI 21 with the sensor facing the floor is recommended. A different mounting position, however, does not affect the sensor's accuracy.

For connecting the CI 21 refer to the connection diagram (also refer to the terminal diagram of your GMA controller). To mount the transmitter, remove the four lock screws and lift the casing top off. Attach the casing by putting two screws through the appropriate borings.



The CI 21 comes equipped with a protective cap over the sensor. Do not remove this cap until just before you apply power to the transmitter.

The protective cap must be replaced whenever power is not applied to the Cl 21.



In order for the CI 21 to detect NH₃, the protective cap must be removed. Once the CI 21 is properly installed, a heater helps to keep the sensor protected from most contaminates in the environment.



If maintenance is being performed in the area, the protective cap must be used.

Failure to follow these instructions may cause damage to the sensor and will void the warranty.



The CI 21 is not certified for use in hazardous areas.

Electrical Connections

An electrician following the applicable regulations should install cable and electrical connections. A shielded 18/3 cable should be used (Belden #8770 or equivalent) for distances less than 500 yards.

NOTE: There is shielding over the EMC cable screwing.




Figure 4


The supply voltage for the transmitter (terminal 4) is 10–30 VDC.

Once the installation is completed, close the casing top and secure it with the screws.

Start-up


Once the CI 21 is installed, the user can put it into operation. Before shipment, the CI 21 is calibrated and tested after a specific warm-up time. When connecting the transmitter to a controller (power supply) you should allow approximately 60 minutes to warm-up the Charge Carrier Injection sensor before checking the sensitivity. Should you notice that a re-calibration is required, allow an extended warm up period of about 4 hours. If you replace the sensor, adjust the jumper of the transmitter according to the specifications on the new sensor (sensor replacement requires adjusting the dead band for the new sensor).

The CI 21 is calibrated and tested before shipment. Certain deviations, however, are possible due to transport, mounting, and ambient conditions. Inside the CI 21 there are two test contacts to check the electrical output signal (fig. 3). The voltage at the test contacts can be measured with a voltmeter. Warm-up time is approximately 60 minutes.



Check Electrical Zero Point

- 1) Use a multimeter to check the output signal by measuring the voltage at the test contacts.
- 2) The reading should be in the range of: 0.198 - 0.200 V
- 3) If there is a deviation from this range it is possible to adjust the value by turning pot 0.




Check Sensitivity (Calibration) / Bump Test

- 1) Use a multimeter to check the output signal by measuring the voltage at the test contacts.
- 2) Attach the calibration adapter to the diffusion inlet of the transmitter. Moisten and wring out the humidifier sponge, install the sponge in the humidifier and attach to the calibration adapter.
- 3) Select a suitable test gas (200, 1,000, or 10,000 ppm NH₃ depending on the detection range) and allow a flow of 0.5 lpm or 1 CFH for 3-5 minutes before reading and making adjustments.
- 4) The reading should be in the 1.00 V +/- 10% range.
- 5) In case of a deviation you can correct the value by turning potentiometer **C**.

When using a GMA controller, signal checks and adjustments may be done at the GMA (please refer to the GMA manual).



Do not use test gas with a balance of nitrogen. Calibration gas must be GfG approved.



Test Contacts

1.00 V


Pot C

0 1 2 3 6 5 4

Check Zero Point (Dead Band)

Adjustments to the zero point (Pot A) should not be done without direct instruction from GfG.

- 1) Use a multimeter to check the output signal by measuring the voltage at the test contacts.
- 2) Attach the calibration adapter to the diffusion inlet of the transmitter. Moisten and wring out the humidifier sponge, install the sponge in the humidifier and attach to the calibration adapter.
- 3) Select a suitable test gas (20 or 30 ppm NH₃, depending on the detection range) and allow a flow of 0.5 lpm or 1 CFH for 3 to 5 minutes before reading and making adjustments.



- 4) The signal should be approximately 0.200 V.
- 5) In case of a deviation it is possible to adjust the value by turning potentiometer **A**.

Pot A is factory set and should not be adjusted without direct instruction or training from GfG.

Troubleshooting			
Symptom	Reason/s	Solution	
0 mA output current	 Incorrect wiring Defective fuse Cut cable/wire DC voltage too high Defective electronics 	 Verify wiring to CI 21 diagram Replace fuse (use approved type only) Reconnect wiring Correct voltage to 10–30 VDC Contact GfG for instructions 	
No reaction to test gas	Improper gas Faulty sensor	Use GfG approved test gasReplace sensor	
Controller shows incorrect reading with test gas applied.	Controller not calibratedCl 21 out of calibration	Calibrate controllerCalibrate CI 21	
High output current with no gas present	Open ground wirePossible interfering background gas	Reconnect ground wireIdentify background gas and eliminate	
Slow sensor response	Incorrect test gasHumidifier sponge too wetHumidifier sponge too drySensor defective	Use GfG approved test gasWring out spongeMoisten spongeReplace sensor	

Service

Service includes maintenance, inspection, and repair of gas warning equipment and function tests of the equipment.

A function test checks the:

- Alarm with gas concentration
- Zero point and sensitivity (calibration)
- Response time
- Gas sampling and gas processing system (if any)
- Alarm signal activation
- Failure alarm

The function check should be done before putting the system into operation, and it should be repeated at least once a year.

For service, please contact:

GfG Instrumentation, Inc. 1194 Oak Valley Drive

Suite 20

Ann Arbor, MI 48108 Toll Free: 800-959-0329 Phone: (734) 769-0573 Fax: (734) 769-1888

E-mail: info@goodforgas.com Website: www.goodforgas.com

Maintenance and Inspection

With maintenance and inspection, you can ensure the proper working status of the gas detection system. The Cl 21 does not require special maintenance. Depending on the ambient conditions, gas detection systems may change their behavior. Therefore, it is important to perform a visual check every day, particularly during the first few days after starting up the device:


- Check pump assembly and gas processing system (if any).
- Check gas supply for blockage (for correct measurement the gas supply to the sensor must not be blocked).
- Sensors are subject to ageing and exhausting. Depending on their type and on their exposure to
 gas they can be used up. A sensitivity check can only be done with a suitable test gas after certain
 periods of time. The sensitivity calibration is an expert calibration and is usually performed by GfG
 service or by authorized technicians.

Spare Parts Part Number Spare sensor/internal Spare sensor/external 737806 737805


If you replace the sensor, adjust the jumper inside the transmitter according to the specifications on the new sensor (sensor replacement requires adjusting the dead band for the new sensor).


Accessories	
	Part Number
Calibration adapter	2000209
Calibration kit without gas (fixed flow regulator, tubing, inline humidifier and adapter)	7721-200
Calibration gas (200 ppm NH₃ balance air)	7803-002
Calibration gas (1,000 ppm NH₃ balance air)	7803-010
Calibration gas (10,000 ppm NH₃ balance air)	7803-100
Calibration kit with gas (10,000 ppm)	7722-100
Calibration kit with gas (1,000 ppm)	7722-010
Calibration kit with gas (200 ppm)	7722-200


Connection Diagram CI 21 (Internal Sensor)




Connection Diagram CI 21 (External Sensor)









Technical Data	
Transmitter CI 21	
Type:	CI 21
Sensor Type:	MK 378-5 (external) or MK 378-6 (internal)
Gas:	Ammonia, NH ₃
Standard range:	200 ppm, 1,000 ppm or 10,000 ppm
Gas Supply:	Diffusion
Detection Principle:	Charge Carrier Injection (CI)
	315 mA, 250 V
Supply Voltage:	
Output Current:	0.2 to 1 mA (4 kΩ) or 4 to 20 mA (200 Ω)
Response Time t _Δ :	
Typical Sensor Life:	> 2 years
Climate Conditions	
Temperature:	-40 to 131°F (-40 to +55°C)
	0 to 99% r. h., non-condensing
Pressure:	920 to 1,100 hPa
Casing	
	Aluminum
Protection:	
	½ inch conduit
Transmitter Cable:	3 core shielded cable AWG 18 for runs less than 500 yards
	Belden #8770 or equivalent
	3.9x3.9x2.2 inches (100x100x57 mm) (HxWxD)
	13 oz. (370 g)
Approvals	
CSA:	c-CSA-us
	CSA C22.2 No. 205-M1983
UL:	UL 916

Warranty

GfG Instrumentation warrants our products to be free from defects in material and workmanship when used for their intended purpose, and agrees to remedy any such defect or to furnish a new part (at the option of GfG Instrumentation) in exchange for any part of any product that we manufacture that under normal use is found to be defective; provided that the product is returned, by the purchaser, to GfG's factory, intact, for our examination, with all transportation costs prepaid, and provided that such examination reveals, in our judgment, that it is defective.

This warranty does not extend to any products that have been subjected to misuse, neglect, accident, or unauthorized modifications; nor does it extend to products used contrary to the instructions furnished by us or to products that have been repaired or altered outside of our factory. No agent or reseller of GfG Instrumentation may alter the above statements.

This page is intentionally left blank.

GfG Instrumentation, Inc.

1194 Oak Valley Dr. Suite 20 Ann Arbor, MI 48108 USA

US/Canada: (800) 959-0329 US/Canada Fax: (734) 769-1888 International: +1 734 769 0573 International Fax: +1 734 769 1888

E-mail: info@goodforgas.com

Website: www.goodforgas.com/us-en